

Model name	Hardware	TSC1 Digital Camera
	Software	RecOrb Imaging Software
System requirements & compatibility	Resolution	4024 x 3036 pixels
	Effective pixels	12.2 MP
	Frame rate	31 fps
	Filter	Yellow filter IN/OUT
	Dynamic Range	70.3 dB (typ.)*
	Compatible models	700GL, 700GL NSW, 700GL LPH, 4ZL, 2ZL
System Requirements	OS	Windows 11
	RAM	≥ 16 GB
	Interface	USB 3.0
	CPU	Intel Core i5 ≥ 3.5 GHz
	Monitor Resolution	Full HD (1920 x 1080)
	Video Board	Graphics memory ≥ 1 GB
	Storage	SSD ≥ 512 GB

*Published value for a monochrome model using the same sensor (EMVA 1288 data for the color model are not publicly available).